



01

## Sterilization



**Product Name:** Disposable protective clothing for medical use

**Product Description:** the product is sterilized by irradiation, and the sterilization period is one month. It can be used by hospital medical staff and also can be used in ICU.

**Product executive standard:** GB19082-2009



02

## Non sterilization



**Product Name:** Disposable protective clothing for medical use

**Product Description:** This product has not been irradiated and sterilized, but it can be used by hospital medical staff, not for ICU.

**Product executive standard:** GB19082-2009

03 **Detailed technology** >>>

## Protective clothing without feet



**Strong sealed brim elastane tape**



**Fly-prssure-sensitive adhesive tape**



**Pressure adhesive tape waist move comfortable**



**Sealed  
pressure  
adhesive  
tape armpit**



**Sealed  
pressure  
adhesive  
tape rise  
and crotch**

## Protective clothing with feet





**X 50**

1 pcs/bag

|                          |                 |
|--------------------------|-----------------|
| <b>Carton Size (cm)</b>  | <b>70*45*45</b> |
| <b>Qty/ctn (box/bag)</b> | <b>50</b>       |
| <b>Gross Weight(kg)</b>  | <b>17.8</b>     |
| <b>Net Weight(kg)</b>    | <b>16.1</b>     |

Dimension table:

## Connected style size

Unit: CM

| Size      | Length | Bust | S.L. | C.W. | Bottom |
|-----------|--------|------|------|------|--------|
| 170       | 173    | 130  | 90   | 18   | 24     |
| 175       | 178    | 135  | 93   | 18   | 24     |
| 180       | 181    | 140  | 96   | 18   | 24     |
| 185       | 188    | 145  | 99   | 18   | 24     |
| Deviation | ±2     | ±2   | ±2   | ±2   | ±2     |