

**PALM FATTY ACID DISTILLATE****Application:**

- Manufacturing of candles, cosmetics, toiletries and pharmaceuticals products.
- Raw material in Vitamin E manufacture as well animal feed manufacturing.

Packing: 185Kg Steel Drum**CAS No: 8002-75-3****Specifications:**

Parameter	Method	Specification
Free Fatty Acid [F.F.A] As Palmitic	AOCS Ca 5a-40	70.0 % Min.
Saponification Value [mgKOH/g]	AOCS Cd 3-25	200.0 - 215.0
Saponifiable matter		95.0 %Min.[Basis 97%]
Total Fatty Matter \ T.F.M 1		95.0 % Min.
Moisture and Impurities	AOCS Ca 2c-25	1.0 % Max
Acid Value	AOCS Cd 3-25	200.0 - 205.0
Iodine Value [mlg/100g%]	AOCS Cd Id-92	55.0 Max.
Titer Deg Â°C		40.0 Min
Colour Lovibond 5 1/4" Max.		4.0 Y 6.0 R
Unsap Matter	AOCS Ca 6a-40	1.5 Max.
Peroxide Value (MEQ/KG)		0.5 Max.
Melting Point		53.0 Max.